



Case Studies

Nancy Mullis, MT(ASCP)SBB

American Red Cross IRL

SEABB 2011

Case 1 - A little bit of everything

- Patient LM, 64 y/o female
- Diagnosis: Not provided
- Referred for antibody ID; history of WAA, anti-K, -Fy^a, -Jk^b.
- Transfusion history:
 - 2 units rbc's 3 months prior to referral
 - 2 units rbc's 1 month ago
 - All units group A, Rh negative

Case 1: Initial results

● ABO/Rh

Anti-A	Anti-B	Anti-D	RhC	A ₁ cells	B cells
4+	0	0	0	0	4+

● DAT

Polyspecific AHG		Anti-IgG		Anti-C3				Saline
IS	Mi	IS	Mi	IS	Mi	RT	Mi	IS/Mi
1+		1+		0	0	0	0✓	0/0

Antibody screen



	IS	LISS 37	LISS IAT		PeG IAT
SI	2+	4+	4+		4+
SII	2+	4+	3+		4+
SIII	1+	2+	2+		3+
AC	1+	0	1+		2+



Phenotype

- We need autologous cells
 - For phenotype to determine what antibodies this patient can produce
 - To investigate if autoantibody is present in the plasma and eluate
- Patient has been transfused within the last 3 months
 - Pre-trx sample not available



Microhematocrit Cell Separation

- Separation of young autologous red cells (i.e. reticulocytes) from transfused RBCs by simple centrifugation in microhematocrit tubes
- Reticulocytes have a lower specific gravity and therefore are concentrated at the top of the microhematocrit tube
- Best if samples collected 3 or more days after transfusion
- Other method: Hypotonic Saline Wash (used for patients with sickle cell disease)

Phenotyping results



- Cell separation performed
 - Used **R**eticulocyte **E**nriched **P**ortion (REP) for testing
- Phenotype:
C-E-c+e+; K-; Fy(a-b+); Jk(a+b-); S-s+;
M+N-; P₁-; Le(a-b-)
- Patient can produce alloantibodies to:
 - **D, C, E, K, Fy^a, Jk^b, S, N, P₁, Le^a, Le^b**
 - History of previous anti-K, -Fy^a, Jk^b

What to do next?



- Is the autoantibody reported by the referring laboratory currently present in the plasma?
- What is reactive at IS and 37° phases?
Could the autoantibody have a broad thermal range, or are these alloantibodies?



Our formula for distinguishing autoantibody from alloantibody

Test DAT negative autologous cells and phenotypically similar cells with patient plasma and/or eluate:

Phen sim pos, auto neg = Probable alloantibody to HIA

Phen sim pos, auto pos = Probable autoantibody

Phen sim neg, auto neg = Probable multiple alloantibodies

Is Warm autoantibody there?

	D	C	E	K	Fy ^a	Fy ^b	Jk ^a	Jk ^b	S	s		IS	LISS 37	LISS IAT	PeG IAT
P S	0	0	0	0	0	+	+	0	0	+		0	0	1+	2+
R E P	0	0	0	0	0	+	+	0	0	+		0	0	1+	2+
C T L	(REP tested with 6% albumin as control)											0	0	0✓	0✓

Antibody screen

	IS	LISS 37	LISS IAT		PeG IAT
SI	2+	4+	4+		4+
SII	2+	4+	3+		4+
SIII	1+	2+	2+		3+
AC	1+	0	1+		2+

Summary so far...




- What we have confirmed in the plasma:
 - Warm-reactive autoantibody reactive at IAT phase in PeG and LISS testing
- What we suspect:
 - Anti-D, -C, -E, -K reactive at IS and 37° phases
- What we have left to do:
 - Confirm -D, -C, -E, -K at IS and 37° (need 2 reactive cells for each specificity)
 - Perform allogeneic adsorption to confirm and r/o common alloantibodies at IAT phase.
 - Perform elution; confirm WAA in eluate

Eluate results

- Presence of autoantibody confirmed in the eluate

	Eluate	Last wash
SI	2+	0✓
SII	2+	0✓
SIII	2+	0✓
REP	2+	0✓
P.S.	2+	0✓



Allogeneic adsorption

- Goal is to remove the autoantibody reactivity by adsorption, while leaving any significant alloantibodies in the plasma.
- Adsorbing cells should lack the antigens that the patient can produce antibody to, so the alloantibody(ies) are not adsorbed.
- We chose to adsorb on D-C-E-K-Jk(b-) red cells after papain-treatment (which destroys MNS and Fy antigens)

Antibody Identification - Selected Cell Panel

Cell #	Rh-Hr					Kell		Duffy		Kidd		Lewis		P	MNSs				Tests with adsorbed plasma (x1)	
	D	C	E	c	e	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s		PeG IAT
1	+	0	0	+	+	0	+	0	0	+	0	0	+	+	+	+	0	+		3+
2	0	0	0	+	+	0	+	0	+	+	0	0	+	+	+	0	0	+		0✓
3	0	+	0	+	+	0	+	0	+	+	0	0	+	+	+	+	0	+		3+
4	0	0	+	+	0	0	+	0	+	+	0	0	+	+	+	+	0	+		2+
5	0	0	0	+	+	0	+	+	0	+	0	0	0	+	0	+	0	+		2+
6	0	0	0	+	+	0	+	0	+	0	+	+	0	+	0	+	0	+		1+
7	0	0	0	+	+	0	+	0	+	+	+	+	0	0	+	+	+	0		1+w
8	0	0	0	+	+	+	0	0	+	+	0	0	+	+	+	+	+	0		2+
9	0	0	0	+	+	0	+	0	0	+	0	0	+	+	0	+	+	0		0✓
10	+	0	0	+	+	0	+	0	+	+	0	0	+	+	0	+	0	+		3+
11	0	+	0	+	+	0	+	0	+	+	0	+	0	0	+	0	0	+		3+
12	0	0	0	+	+	0	+	0	+	+	0	+	0	0	0	+	+	0		0✓

Testing done to finish up...



- Additional testing was completed with neat plasma that confirmed presence of anti-D, -C, -E, and -K at IS and 37^o phases of testing.
- Additional testing was completed with adsorbed plasma that confirmed the presence of anti-D, -C, -E, -K, -Fy^a, and -Jk^b at the IAT phase of testing.
- The presence of anti-S and anti-N was excluded.



Hmmm, something doesn't make sense....

- The patient was only transfused with Rh negative rbc's
 - Confirmed by referring facility
 - Patient received no platelets or RhIG products
- Where did the anti-D come from???

The G antigen



- G antigen is due to an amino acid that is present on Rh proteins that express the C or D antigen
- The G antigen is rarely expressed on D-C- red cells (r^G and r''^G)
- There are also rare examples of D+ cells that do not express G

Antibody Identification - Selected Cell Panel

Cell #	Rh-Hr					Kell		Duffy		Kidd		Lewis		P	Other typings	Neat plasma	
	D	C	E	c	e	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁		IS	LISS 37
1	0	0	0	+	+	0	+	0	+	+	+	0	+	0	D-C-, G+	2+	3+
2	0	0	0	+	+	0	+	+	+	+	+	0	+	+	D-C-, G+	2+	3+
3	0	0	0	+	+	0	+	+	0	+	0	0	0	+	(G-)	0	0
4	0	0	0	+	+	0	+	0	+	0	+	0	+	+	(G-)	0	0

- Presence of Anti-G confirmed at IS and 37° phases of reactivity
- The presence of anti-C and -D were not specifically excluded
 - Would require use of adsorption studies
 - Academic interest only, except in cases to evaluate need for RhIG

Conclusions



- Antibodies present in plasma:
 - Anti-G reactive at IS, LISS 37°
 - Anti-K, -E reactive at IS, LISS 37°, PeG-IAT
 - Anti-Fy^a, -Jk^b reactive at PeG-IAT
 - Warm-reactive autoantibody reactive by PeG-IAT and LISS-IAT
 - Not ruled out: Anti-D and -C
- DAT positive (IgG); Eluate: Panagglutinin
- Trx. Recommendation: ABO compatible rbc negative for D, C, E, K, Fy^a, and Jk^b, crossmatch compatible with neat or adsorbed plasma.

Case 2 - Do you like to gamble?

- Patient MS, 21 y/o male
- Diagnosis: Sickle cell crisis
- Referred for STAT Antibody identification with 4.8 hemoglobin
- Transfusion history
 - 2 units rbc's transfused 8 days prior to sample
 - No other history given

Case 1: Initial results

● ABO/Rh

Anti-A	Anti-B	Anti-D	RhC	A ₁ cells	B cells
0	0	4+	0	4+	4+

● DAT

Polyspecific AHG		Anti-IgG		Anti-C3				Saline
IS	Mi	IS	Mi	IS	Mi	RT	Mi	IS/Mi
1+		1+		0	0	1+		0/0

Antibody screen and Eluate

	IS	LISS 37	LISS IAT	PeG IAT		Eluate	Last wash
SI	0	0	0✓	2+		3+	0✓
SII	0	0	0✓	2+		3+	0✓
SIII	0	0	0✓	2+		3+	0✓
AC	0	0	1+	1+		NT	NT

LISS panel

Rh - Hr							Kell				Duffy		Kidd		Lewis		P	MN				Lutheran		Xg	PATIENT'S SERUM TEST RES TEST METHODS						
D	C	c	E	e	V	C ^w	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a					
+	+	0	0	+	0	+	+	+	0	+	0	+	+	0	+	+	0	+	+	+	+	0	0	+	0	+	0	1	0	0	✓
+	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	0	+	0	+	0	+	+	2	0	0	✓
+	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	0	+	0	+	+	3	0	0	✓	
+	0	+	0	+	+	0	0	+	0	+	0	+	0	0	+	+	0	0	+	0	+	0	0	0	+	+	4	0	0	✓	
0	+	+	0	+	0	0	+	+	0	+	0	+	+	+	+	+	0	0	0	+	+	+	0	0	+	+	5	0	0	✓	
0	0	+	+	+	0	0	0	+	0	+	0	+	0	0	+	0	0	+	+	+	0	+	+	+	+	6	0	0	✓		
0	0	+	0	+	0	0	+	+	0	+	0	+	0	+	+	+	+	0	0	+	0	+	+	0	+	+	7	0	0	✓	
0	0	+	0	+	0	0	0	+	0	+	0	+	+	0	+	+	+	0	0	+	+	0	+	0	+	0	8	0	0	✓	
0	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+	0	0	0	+	0	+	+	0	0	+	+	9	0	0	✓	
+	+	0	0	+	0	0	0	+	0	+	0	+	0	+	+	0	+	0	0	+	0	+	0	0	+	+	10	0	0	✓	
+	+	0	+	+	0	0	0	+	0	+	0	+	+	+	+	+	0	0	+	+	0	+	0	0	+	+	TC	0	0	✓	

LKSS
37146 u

Can we interpret these results?

- Only one small EDTA sample was sent from the referring hospital; only enough plasma left for crossmatching.
- No red cells remaining for cell separation and phenotype (used to prepare eluate).
- This looks like a warm autoantibody; alloantibodies are ruled out by LISS testing. Should we send units because this is a STAT, or ask for more sample to complete the investigation?

More work to do...

- Received additional sample
- Need autologous cells and phenotype
 - Patient recently transfused
- Confirm presence of autoantibody

Phen sim pos, auto neg = Probable alloantibody to HIA

Phen sim pos, auto pos = Probable autoantibody

Phen sim neg, auto neg = Probable multiple alloantibodies

Phenotyping results



- Cell separation performed
 - Used **Hypotonic Saline Wash method**
 - 0.3% saline lyses normal donor rbc's but leaves HbSS cells intact
- Phenotype:
C-E-c+e+; K-; Fy(a+); Jk(b+); S-s-; U-
Fy^b and Jk^a typings inconclusive (weaker than expected reactions)
- Patient can produce alloantibodies to:
 - C, E, K, [Fy^b, Jk^a], S/s/U

Is Warm autoantibody there?

	D	C	E	K	Fy ^a	Fy ^b	Jk ^a	Jk ^b	S	s	IS	PeG IAT	Eluate
U+	+	0	0	0	0	0	0	+	0	+	0	2+	3+
U- PS	+	0	0	0	0	0	0	+	0	0	0	0✓	0✓
REP	+	0	0	0	0	?	?	+	0	0	0	0✓	0✓
CTL	(REP tested with 6% albumin as control)										0	0✓	0✓

Summary



- Additional U- cells were tested with plasma and eluate to confirm anti-U and rule out other common allos by PeG-IAT.
- A delayed transfusion reaction (with antibody to HIA or to multiple antigens) can easily be confused serologically and clinically with WAIHA
 - Positive DAT; may appear mixed-field, but MF is hard to detect if the reactivity is weak
 - All cells reactive with plasma and eluate, including the AC if the patient's rbc sample includes donor rbc
 - Signs of hemolysis - may include unexplained drop in Hct/Hgb, jaundice, hemoglobinuria, etc.
- Transfusion history is important!

Happy Spring!

